

## Sample Results Summary Sheet

Please return this form to the Curator for each allocated Sample

Sample ID: RA-QD02-0049

PI: Hiroshi Naraoka

**Type and date of analysis performed:**

Amino acid analysis for organic solvent extract by two-dimensional high performance liquid chromatography with highly sensitive fluorescence detection.

**Elements or phases identified:** (Mg, Si, olivine, pyroxene, aromatic carbon, etc.)

Not determined

**Contaminant phases identified:** (Al, SUS, carbon particles, etc.)

Not determined

**Sample handling:** (e.g. exposed in atmosphere, embedded in resin, polished, sliced by FIB or UMT)

Extracted with dichloromethane/methanol (1/1, ~0.2ml x3) on a clean bench under atmosphere

**State of sample pre-analysis:** (e.g. N2 hold, atmosphere, resin embedded, polished section, UTS) (please describe treatments and/or modifications for the sample you have done before your analysis)

In a diamond holder after Raman and Infra-Red spectroscopy measurement

**State of sample post-analysis:**

(N2 hold in sample holder, atmosphere, resin embedded, polished section, UTS)

(partially damaged by electron beam, spotted by Ga beam, neutron activation)

(consumed by laser ablation)

(unexpected breakup, into # pieces)

(Lost : reason)

Resin embedded and stored under N2 for Synchrotron-Radiation tomography

**Analysis data Notes:** (summary of the attached analysis data and/or images)




Fig. 4. Amino acid chromatograms of glycine (Gly), alanine (L,D-Ala),  $\alpha$ -aminoisobutyric acid (AIB) and isovaline (L,D-Isoval). a) standard, b) RA-QD02-0049, c) RA-QD02-0033, d) blank, and e) Y-791191. Samples of RA-QD02-0049, RA-QD02-0033, blank and Y-791191 were injected by 40  $\mu$ l of whole solution (100  $\mu$ l).




Fig. 4. (continued).